

## N SA DICA Nº A NEW SANDAY VILLE GE

AN VIT AN MITOMATERE ERRESEVER SHVATE COMES:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

November 05, 2004

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.

> APPLICATION NUMBER: 60/510,237 FILING DATE: October 10, 2003

RELATED PCT APPLICATION NUMBER: PCT/US04/33888

Certified by

Jon W Dudas

Acting Under Secretary of Commerce for Intellectual Property and Acting Director of the U.S. Patent and Trademark Office



PTO/SB/16 (08-03) Approved for use through 07/31/2006. OMB 0651-0032 U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

PROVISIONAL APPLICATION FOR PATENT COVER SHEET This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c). Express Mall Label No. EV329487945US

F		INVENTOR	(S)				]
Given Name (first and middle [if any])		Family Name or Surname		(City ar		Residence State or Foreign Country)	ļ
Dongqing		Chen		Port Jeffer	son Static	on, New York	1
Sarang				Selden, No			4
Additional inventors are be			separately numb		ttached he	ereto	4
		LE OF THE INVENTION (			Calcal Ca		ᅥ
Method and System for Direct all correspondence	Preparing 3-Dime	nsional Volumetric Image	s of a Selected	1 Organ for v	intual En	doscopy	ᅦ
Customer Number:							١
Customer Number.							ı
OR							╛
Firm or Individual Name	F. Chau & Associate	es, LLP					_
Address	1900 Hempstead Tu	ırnpike					긕
Address	Suite 501				7!-	T	ᅴ
City	East Meadow		State	NY	Zip	11554	4
Country	U.S.		Telephone	46-4	Гах	516-357-0092	닉
	ENCLO	SED APPLICATION PAR	KTS (check all	tnat apply)		<del></del>	ㅓ
Specification Numb	er of Pages 27	<del></del>		CD(s), Numbe			
Drawing(s) Number	of Sheets		V	Other (specify)	Form P	TO-2038	
Application Date Sh	eet. See 37 CFR 1.7	<b>'</b> 6					
METHOD OF PAYMENT	OF FILING FEES F	OR THIS PROVISIONAL API	PLICATION FOR	PATENT			
Applicant claims sr	nall entity status. See	e 37 CFR 1.27.		FILING FEE			
	order is enclosed to				Amou	unt (\$)	
					] .		
The Director is her	The Director is herby authorized to charge filing fees or credit any overpayment to Deposit Account Number:				(	80.00	
Payment by credit card. Form PTO-2038 is attached.							
The invention was made by an agency of the United States Government or under a contract with an agency of the							
United States Government.							
V No.							
Yes, the name of the	ne U.S. Government	agency and the Government	contract number	r are:			
		[Page 1 c	of 2]	Date October	10, 2003		
Respectfully submitted						13 584	
SIGNATURE	SIGNATURE TO STATE TO			(if ennmnriate)	EGISTRATION NO. 43,584 appropriate)		
TYPED or PRINTED NAME Frank V. DeRosa				Docket Number: 8095-7			

TELEPHONE (516) 357-0091

USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT

USE ONLY FUR FILING A PROVISIONAL APPLICATION FOR PATENT

This collection of information is required by 37 CFR 1.51. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 7 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden, should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Mail Stop Provisional Application, Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

#### PROVISIONAL APPLICATION COVER SHEET **Additional Page**

PTO/SB/16 (08-03)

Approved for use through 07/31/2006. OMB 0651-0032

U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

Docket Number 8095-7 INVENTOR(S)/APPLICANT(S) Residence (City and either State or Foreign Country) Given Name (first and middle [if any] Family or Surname Sunnyvale, California Kreeger Kevin Greenlawn, New York Wax Mark Plainview, New York Kaufman Агіе Liang Stony Brook, New York Zhengrong

Number

[Page 2 of 2]

## METHOD AND SYSTEM FOR PREPARING 3-DIMENSIONAL VOLUMETRIC IMAGES OF A SELECTED ORGAN FOR VIRTUAL ENDOSCOPY

Inventors: Dongqing Chen, Kevin Kreeger, Sarang Lakare, Mark Wax, Arie Kaufman, and Zhengrong Liang

Assignee: Viatronix Incorporated, Stony Brook, New York

#### **References Cited**

The following references are fully incorporated herein by reference.

#### U.S. Patent Documents

5,782,762	7/21/1998	Vining et al.
5,971,767	10/26/1999	Kaufman et al.
6,331,116 B1	12/18/2001	Kaufman et al.
6,343,936 B1	02/05/2002	Kaufman et al.
6,477,401 B1	11/05/2002	Johnson et al.
6,514,082 B2	02/04/2003	Kaufman et al.

#### Other Publications


- 1. Wax et al.: "Electronic colon cleansing for virtual colonoscopy," 1998, The 1st International Symposium for Virtual Colonoscopy, Boston, MA, p. 94.
- 2. Chen et al.: "A novel approach to extract colon lumen from CT images for virtual colonoscopy," 2000, IEEE Trans. on Med. Imaging, vol. 19, p. 1220-1226.
- 3. Chen et al.: "A multi-scan MRI-based virtual cystoscopy," 2000, Physiology and Function from Multidimensional Images, Chin-Tu Chen, Anne V. Clough, Editors, Proceedings of SPIE vol. 3978, p. 146-152.
- 4. Lakare et al.: "Automated pre-navigation processing for virtual colonoscopy," 2000, The 2<sup>nd</sup> International Symposium on Virtual Colonoscopy, Boston, MA.
- 5. Wax et al.: "Optimizing bowel preparation for virtual colonoscopy electronic cleansing," 2001, The 87<sup>th</sup> Scientific Assembly and Annual Meeting of RSNA, Chicago, IL.
- Callstrom et al.: "CT colonography without cathartic preparation: Feasibility study," 2001, Radiology, vol. 219, p. 693-698.
- Zalis et al.: "Digital subtraction bowel cleansing in CT colonoscopy," 2001, Am J. Roentgenol, vol. 176, p. 646-648.
- 8. Chen et al: "Laxative free virtual colonoscopy: Feasibility study," 2002, The 3<sup>rd</sup> International Symposium on Virtual Colonoscopy, Boston, MA.
- 9. Wax et al: "Optimizing laxative free virtual colonoscopy," 2002, The 88th Scientific Assembly and Annual Meeting of RSNA, Chicago, IL.


#### Abstract

A method for preparing 3-dimensional volumetric images of a selected organ for virtual endoscopy is provided. The patient is administrated either with an amount of contrast agent, or specific diet, or combination of both which will tag the residue in the lumen of the organ and allow the residue to be distinguished from surrounding tissues in the medical images. The patient's organ is imaged by computed tomography (CT), or magnetic resonance imaging (MRI), or other available modalities after the management of diet and contrast agent. The acquired 3-dimensional volumetric images are then processed with tissue and organ segmentation and feature analysis focusing on the region of selected organ. The extracted feature data of clinical significant is transformed and fused with the original volumetric images for either removing residue in the lumen or facilitating the contrast between normal and abnormal

tissues. The fused volumetric images are displayed and rendered for physician inspection of abnormalities. Applications to laxative free CT virtual colonoscopy are provided.

#### **Figures**








Figure 3.

FOOD GROUPS	FOOD ALLOWED	FOODS TO AVOID
Breads	White refined breads, rolls biscuits, muffins, crackers, pancakes and waffles	Whole grain flower products of any type and baked goods made with bran, nuts, seeds, coconut, fruits, bagels, cornbread or graham crackers
Cereals	Refined cooked cereals, including cream of wheat, and farina, puffed wheat, puffed rice and rice krispies	Oatmeal, any whole grained cereal bran or granola and any containing nuts, seeds, coconut or dried fruit.
Desserts	Plain cakes and cookies, water ices (Marinos), plain low-fat yogurt, Jello, custard, grape or apple jelly, plain hard candy, marshmallows and lite ice cream without nuts or chocolate	Any desserts made with whole grain flour, bran, seeds, coconut, dried fruit, yogurts with fruit skins or seeds or nuts, sherbets and popcorn. No chocolate.
Fruits	Most canned or cooked fruits, applesauce and ripe banana.	Dried fruits, all berries, most raw fruit, except banana.
Potato and Potato Substitute	Cooked white potato without the skin, white rice, white pasta and egg noodles	All others, including whole-wheat pasta, noodles, vegetable pastas, and sweet potato.
Vegetables	Cooked fresh, frozen, or canned carrots, beets, asparagus tips, French style string beans and pureed squash, spinach.	All others including raw and deep-fried vegetables, broccoli, cauliflower cabbage, spinach, sauerkraut, winter peas, corn, and any other vegetables with seed.
Fats	Margarine, salad oil, lite salad dressings, lite mayonnaise, and plain gravy.	Butter, any fat containing whole grain flour, bran, seeds, nuts, Coconut, or dried fruit.
Meats and Meat Substitutes	Ground and well-cooked white meat chicken and turkey, with skin removed, Fish, shellfish, eggs, and low-fat cheese.	Red meat, BBQ or pickled meat, any made with whole grain ingredients, seeds or nuts, dried beans, peas, lentils, legumes, peanut butter and whole milk Cheese.
Soups	Bouillon, broth, low-fat cream made with allowed vegetables, noodles, rice, or refined white flour.	All others.
Beverages	Decaffeinated liquids of all kinds, Caffeinated beverages limited to 2-3 (10 oz.) cups per day, low- fat milk, and strained fruit juices.	Espresso, frappucino, cappuccino, whole milk, fruit/vegetable juices containing pulp, prune juice, and all alcoholic beverages.

Figure 4.

Foods in the provided Meal Kit
Vanilla Nutritional Shakes
Lemon drinks
Chocolate flavored energy bars
Cinnamon apple sauce
Potato poppers
Instant noodle soup
Stroganoff

Figure 5.

Food Groups	Food Allowed	Foods to avoid
Beverages	Apple, Welch's white grape juices,	All other.
	Gatorade (no red, orange or purple) Decaffeinated weak tea.	
Soups	Clear Broths or Bouillon or consommé.	All others
Desserts	Lemon or lime flavored Jell-O, water ices and frozen ice pops	All others

Figure 6.




Figure 7.

	Num. of Cases	Num of	good tagging cases	Num. o	of poor tagging cases
BP 1	30	30	100%	0	0%
BP 2	35	31	88.57%	4	11.43%
BP 3	14	7	50%	7	50%

Figure 8.




Figure 9.




Figure 10.




Figure 11.




Figure 12.




Figure 13.





Figure 14



	Nun.	<b>Electronic Cleansing Quality</b>			
		Adequate Inad		lequate	
Diet 1	20	20	100%	0	0%
Diet 2	26	20	77%	6	23%
Diet 3	14	8	57%	6	43%
Total	60	48	80%	12	20%

Figure 16.




Figure 12. This picture demonstrates an example of polyp rendering in different color in the virtual colonoscopy endoluminal view.

### METHOD AND SYSTEM FOR PREPARING 3-DIMENSIONAL VOLUMETRIC IMAGES OF A SELECTED ORGAN FOR VIRTUAL ENDOSCOPY

#### FIELD OF THE INVENTION

The presented invention relates generally to a method or system for preparing 3-dimensional (3-D) volumetric images of a selected organ for virtual endoscopy. The virtual endoscopy is generally for early cancer detection related to organ with hollow lumen. One such application is for generating CT virtual colonoscopy images that does not require applying neither oral laxative nor suppository to the patient bowel preparation.

#### **BACKGROUND OF THE INVENTION**

Virtual endoscopy has received increasing attention as a minimum-invasive and patient comfortable method for evaluating hollow organs. The fast increasing of virtual colonoscopy practice is an evidence that the virtual endoscopy approach is getting better acceptance. The basic idea of this technology is to fill the selected organ with materials that can be distinguished from the wall of the organ in the medical images. The acquired medical images are further processed for constructing either 2-D or 3-D model of the organ. The model is displayed in different rendering modes for physician's inspection of abnormalities (US Patents 5,782,762 and 5,971,767).

The hollow organ may not show any cavity in the medical images if it is empty. For example, the empty colon and bladder usually collapse and show no cavity in CT images. Even if the lumen of organ is filled with some material, the lumen may still not be able to be distinguished from surrounding tissue. For example, the bladder that is filled with urine may not be distinguished well from surrounding tissue in CT images. Furthermore, there may be more than one kind of materials inside the lumen of the organ. For example, there are usually air, stool, and fluid in the human colon. That results more difficulties to detect abnormality related to the organ. Hence, for virtual endoscopy applications, it is necessary to manage either a certain amount of contrast agent (oral intake, suppository, or intravenous injection), or diet, or both of them to the patient for achieving optimal imaging results.

The procedure of the administration of both contrast agent and diet is called organ preparation. The organ preparation may take several days or several minutes. The main concern is to achieve optimal tissue contrast in the images. Another concern about the organ preparation is non-invasive and patient comfortable. For example, there have been evidences showed that the major complaint to the optical colonoscopy is the bowel cleansing rather than the catherter procedure itself. Although the virtual colonoscopy is free from perforation, the physical cleansing bowel preparation procedure is one of the major barrier to allow virtual colonoscopy to be a screening tool for colorectal cancer. It is highly desirable for an easy handle, non-invasive, and patient comfortable bowel preparation for virtual colonoscopy.

There are usually several hundreds of axial images for a single virtual endoscopy scan. It is a time-consuming work to review all those axial images. The virtual endoscopy usually provides a 3-D view of the inner surface of the organ. It is called endoluminal view. The physician can either change viewing and lighting source and direction or

automatically navigate inside the organ. The 3-D display organizes and shows the anatomy information in the slice images in a more efficient way. There have already been a lot of approaches for displaying and rendering the virtual endoscopy volumetric images (eg. US Patent, 5,971,767, US Patent 6,331,116, US Patent 6,477,401, and US Patent 6,514,082). Those approaches basically focus on the mode of navigation inside the lumen and the speed of the rendering. They are efficient for shape description. If there is an irregular bump shooting out from the inner surface of the organ (like the polyp in the air-distended colon), it is easily detected in the virtual endoscopy. However, there are abnormalities happening right in the wall of the organ and resulting only mild shape deformation on the inner surface, especially for lesion at the early stage (like the bladder lesion and cancer). As in the MRI virtual cystoscopy images, the lesion that invading the wall of the organ demonstrate different texture comparing to the normal wall. If the texture information may be overlaid on the endoluminal view, it will greatly facilitate the early detection of lesion and cancer. If the texture information may be fused in the original volumetric images in an appropriate way, it is possible that the lesion region will be enhanced in the endoluminal view. Hence, the method for fusion volumetric images and feature data is highly desirable. This fusion technology should also be very helpful for displaying computer assisted lesion detection result in its real anatomical context.

#### SUMMARY OF THE INVENTION

.....

#### BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1 is a flowchart representing a method in accordance with the present invention of preparing volumetric images of a selected organ for virtual endoscopy.

Figure 2 is a flowchart representing a method and system of laxative free virtual colonoscopy.

Figure 3 is a chart illustrating in detail the administration schedule and doses of the contrast agent for laxative free virtual colonoscopy.

Figure 4 is a chart describing the suggested lower residue diet. It tells the patient what kind of food should be avoided and what kind of food can be taken during the bowel preparation period.

Figure 5 is a chart describing the foods that the Meal Kit provides.

Figure 6 is a chart describing the foods that can be taken for the fluid diet.

Figure 7 is a comparison chart of those three diets for laxative free virtual colonscopy.

Figure 8 is a chart to list tagging quality of those three bowel preparations for laxative free virtual colonoscopy.

Figure 9 demonstrates electronic cleansing results on virtual colonoscopy images. The tagged residue (bright region) that is shown in the left panel in the colon lumen is removed in the right panel (becoming dark as air region). The layer between tagged residue and surrounding tissue in the left image is replaced by the layer between air and tissue in the right image. The region in the circle demonstrates a tiny tagged stool attaching on the colon wall. The tagged residue puddles that are arrowed show the non-uniform tagging.

Figure 10 is a flowchart describing steps of the electronic cleansing method.

Figure 11 shows a CT image that the colon wall has the similar intensity value as that of tagged residue. The arrowed 719 HU shading region is a thin haustral fold. The arrowed brighter region of 823 HU is a puddle of tagged fluid in the colon.

Figure 12 a chart that describes the 5 directional gradients for extraction of tissue boundary in CT images.

Figure 13 is a chart showing the neighbor voxels that are used in computation of selected directional gradients for CT images.

Figure 14 demonstrates the maximum gradient image (right) of a CT axial image (left).

Figure 15 demonstrate the result of electronic cleansing in both 2-D and 3-D display. The top left image is the original CT slice image. The thin layer of tagged stool attaching to the ceiling of the colon is encircled. The top right image shows the result of electronic cleansing. The bottom left image is the 3-D view of the colon wall that is covered by the thin layer of tagged stool. The bottom right image is the 3-D view of the colon wall after removal of the thin layer of tagged stool. It is obvious that the small thin haustral folds are better demonstrated in the cleansed 3-D display.

Figure 16 lists the quality evaluation results of the invented electronic cleansing algorithm. The result is based on the radiologist's reviewing of 60 cases from 35 healthy volunteers.

Figure 17 demonstrates result of fusion between the computer assisted poly detection and the original volumetric images. The detected polyp shows different color from that of normal colon wall in the endoluminal view.

#### DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENTS

The present invention generally relates to a method and system, as schematically represented in Fig. 1, for preparing 3-D volumetric images of a selected organ for virtual endoscopy. While the methods and system described in this application can be applied to any object with hollow lumen to be examined, the preferred embodiment is an body organs with hollow lumen such as colons, tracheo-bronchial airways, bladder and the like.

The following is the detail description of the invented technologies for laxative free virtual colonoscopy.

The method and system for laxative free virtual colonoscopy is described schematically in Fig.2. It illustrates major steps to perform a laxative free virtual colonoscopy. The result of polyp detection represents all suspicious regions or areas in the images that are detected automatically by algorithms. Whether the result of polyp detection should be fused with the original volumetric images depends of the decision of physician who uses the virtual colonoscopy. If the result of polyp detection is not available, the physician can still use the virtual colonoscopy system for polyp detection.

The preferred embodiments including design of low residue diets and administration of doses of contrast agents for a laxative free bowel preparation for virtual colonoscopy that provides good residual stool and fluid tagging quality. The preferred embodiment also including intelligent algorithms for electronic cleansing of tagged residual stool and

fluid in the virtual colonoscopy images. The electronic cleansed volumetric images are then displayed for physician inspection.

We use barium sulfate (2.1%, 250 ml) and Gastroview (diatrizoate meglumine and diatrizoate sodium solution, 367 mg/ml, 120 ml) as colon residue contrast agent. The barium sulfate is in banana flavor. The Gastroview can be mixed with soft drink, like soda, if the patient feels more easily to take. The optimal dosage to administrate the contrast agent is the 36 hours 6 doses as described in the Fig. 3. Neither oral laxative nor suppository is given to the patient. There are several kinds of barium and sodium iodine available in the market. They result similar tagging quality with the same dosage. Hence, the invention is not limited to specific brand named contrast agent.

Without applying laxative, the stool and fluid residue may remain in the colon anywhere. The patient is required to follow specific diets within 3 days in order to facilitate uniform residue tagging and to mitigate the stool sticking around the colon wall. The design of diets focuses on easy handling and less interfering of patient's normal daily life. There were three kinds of three-days-diet being designed. The first diet (denoted as Diet 3) was a general lower residue diet. It was very close to the normal diet. If patient follows Diet 3, he/she is only to avoid several kinds of food that possibly result more residue for 3 days prior to imaging. The detail food list for the Diet 3 is shown in the Fig. 4. The second diet (denoted as Diet 2) let the patient following the same food list in day 3 and 2 prior to the day of scan and required the patient eating food only from a provided Meal Kit in the day 1 prior to the day of scan. The foods that provided in the Meal Kit are shown in the Fig. 4. The third diet (denoted as Diet 1) required the patient following a liquid food diet in the day 1 prior to the scan. In the other 2 days, the requirement of the Diet 1 is the same as that of the Diet 3. The foods for the fluid diet are shown in the Fig. 6. The Fig. 7 is a comparison matrix for the selected three kinds of diets. During the 3-days diet for all those 3 kinds of diets, the patient was suggested to take small portion of meal and avoid the following foods: whole grain flour or cereals, dried or raw fruits, raw or deep fried vegetables, tough fibrous meats, caffeinated liquids, nuts and seeds, and yogurt. The bowel preparation procedure (BP) for virtual colonoscopy combines one of the 3-day diets and the administration of 6 doses contrast agent. The bowel preparation procedure associated Diet 1, 2, and 3 are denoted as BP 1, 2, and 3 respectively.

Prior to the scan, approximate 2 to 3 liters room air or CO2 was forced into patient's colon in order to expand it to make the colon easier to scan and exam. The routine CT virtual colonscopy scanning protocol was utilized. An example of a helical CT scanner (single slice) would use 3 to 5 mm collimation, 1:1 to 2:1 pitch, 120 kVp, and 100 mA. The patient was required to hold a breath. The scan usually took 35 to 45 seconds on the helical CT scanner depending on the body size of patient. It generated 300 to 550 images for each scan series. The images were packed to form a volumetric dataset for further processing or reviewing. Usually, the patient was scanned twice in different body positions at a single visit to scanning room. One time the patient laid down with face up (supine position) and another time the patient laid down with face done (prone position).

The effectiveness of the bowel preparations was evaluated on 45 volunteers (average age 36, standard deviation 15, minimum 18, maximum 69, 6 females and 39 males). Each volunteer took at least one of the bowel preparations and underwent CT in one day. Some of them took two or all 3 of those bowel preparations and underwent CT scans in

different days (at least one week later). The acquired images were reviewed slice by slice by an experienced radiologist to evaluate the tagging quality. If intensity of residue voxel that was less than 150 Hounsfield (HU) Units was found in the images, the residue was defined as poorly tagged. If one region of residue that was larger than 5 mm in diameter was found poorly tagged in the volumetric images, the dataset or the case related to the volumetric images was defined as poor tagged. Fig. 8 showed the results of tagging quality on those volunteer datasets. There were 79 datasets in total. Among them, there were 30 from BP 1, 35 from BP 2, and 14 from BP 3. From the Fig. 8 we can conclude that the BP 1 achieved the best tagging quality.

The goal of electronic cleansing is to detect and remove all tagged residue inside colonic lumen and to rebuild the image of boundary between tagged residue and colon wall into the image of boundary between air and colon wall. The Fig. 9 demonstrates a result of electronic cleansing. There are a lot of factors that make the electronic cleansing a tough task for laxative free virtual colonscopy CT images. First, the intensity range of tagged residue is similar to that of bone in CT images. The tagged residue in the rectum and sigmoid may frequently touch the backbone in supine position. That makes the separation of tagged residue from backbone difficult. Second, the residue may also be tagged in a non-uniform way so that a parameter set for cleansing one puddle successful may fail for cleansing another. Third, some of thin haustral fold submerged in the tagged residue may have very similar intensities that other tagged residue has due to partial volume effect of CT imaging. The previous inventions in the US Patent 6,331,116, US Patent 6,343,936, and US Patent 6,514,082 provided electronic cleansing algorithms for virtual colonoscopy CT images. However, it was assume that the residue was tagged uniformly and formed a puddle with flat level due to the gravitation. The technologies there could not cleanse tagged stool attaching around the colon wall that does not form a puddle. Unfortunately, the stool often attaches around the colon wall if no laxative or suppository is applied. In this invention, new methods were presented to achieve better quality of electronic cleansing for any kinds of residue tagging condition and any morphology of residue region.

The presented electronic cleansing method was composite of several steps. Fig. 10 shows a flowchart that representing those steps. The first step of electronic cleansing was the extraction of the colon region. In the laxative free virtual colonoscopy images, the colon region was composite of air, residue, and wall. The thickness of colon wall was within a certain range, e.g. less than 5 mm. If the colon lumen, including both air and residue, was found, the colon wall region could be determined by dilation the region of colon lumen. In the CT images, intensity of bone is similar to that of tagged residue and the lung is also filled of air. In US Patent 6,331,116, US Patent 6,343,936, and US Patent 6,514,082, methods had been provided for detecting and removing the regions of both lung and bone while keeping and labeling the region of colon lumen. In this invention, the detection, cleansing, and boundary reconstruction of the region of tagged residue are the focuses.

When all tagged residue regions were found, the low-level voxel classification algorithm that invented in the US Patent 6,331,116 was applied to all regions of tagged residue. The voxels in those regions were grouped into several clusters. Usually, there were 10 clusters being generated. The voxels with similar intensity properties were in the same group. For example, voxels that located in the center of uniform tagged region were in the same cluster and

voxels that located in the edge of an uniform tagged region were in another cluster that was different from that of the center cluster. However, the classification results basically represented intensity properties of the entire volumetric images. The voxels in the same cluster might represent different clinical meaning due to be in different tagged regions with different tagged conditions. For example, the intensity of the fold that was covered by tagged fluid might have similar intensity of tagged fluid due to partial volume effect. Hence, electronic cleansing algorithm must be able to self-adaptive to variation of tagging condition in different tagged regions.

The basic idea of the new algorithm was based on the following fact. In virtual colonoscopy images, the main concern was to detect the boundary between different tissues or materials. In laxative free virtual colonoscopy images, the concern was to detect and extract boundaries between air, colon wall, and tagged residue. Although the intensities of air, colon wall, and tagged residue were in separate ranges that separate from each other, the boundaries of them in the CT images could not be extracted using the simple thresholding due to non-uniform tagging and partial volume effect in the CT images. For example, the colon wall tissue often has intensity around [-100, 50] HU and the well-tagged residue has intensity around [200, 800] HU. If a thin haustral fold submerges in well-tagged residue puddle, the intensity of the colon wall of the thin fold may greater than 600 HU as the example shown in the Fig. 11. The immediate solution was the edge detector. Since most edge detection technologies were based on the difference of intensity rather than the intensity range itself, they were good candidate for determining the boundaries. However, the conventional edge detection technology could not be applied directly here. First, the volumetric CT images often have anisotropic voxel size. Of course, the images can be re-sampled or interpolated into isotropy voxel size. But that takes extra computation time. Second, the task for electronic cleansing is to extraction of the layer of the boundary rather than line of edge curve. The layer of the boundary represents the partial volume effect, which comes from the physical limitation of CT scanner. Third, it is well known that most edge detection method is sensitive to noise and the noise level of CT images usually is high.

In the current invention, a boundary layer extraction method was provided which was based on the maximal gradient approach. The gradient is defined as the directional first-order derivative. Usually, the integer grid of medical images has different unit lengths along different coordinate directions due to the variation of scanning protocols. When computing the discrete first-order derivative in the digital images, the anisotropic voxel size must be considered. Another concern is the so-called partial volume effect. The tissue boundary usually forms a thick layer with certain range due to the partial volume effect rather than a sharp line of curve. The intensity of voxels in the boundary layer changes from one type to another. The boundary layer should include more information on tissue abnormality or tissue feature rather than that of a sharp boundary line does. The maximal gradient approach is design to meet these purposes. The basic idea of the approach is to selected several discrete directions in the 3-D image grid. These directions should be uniformly covered all directions. Then, compute first-order derivatives along all those directions for each voxels and determine the maximum of the absolute value of those first-order derivatives associated. The maximum value is denoted as the gradient feature value (GFV) of that voxel. By thresholding the GFV, the boundary layer can be determined.

For simplicity, we describe the algorithm with a CT image that has voxel size around X:Y:Z=0.7:0.7:1.0 mm. For each voxel in the image, we compute five directional derivatives. Fig. 12 shows the five directions. The X, Y, and Z are the orthogonal axial directions of the image grid coordinate system. The other two directions are the diagonal directions on the X-Y plane. In this invention, the directional derivative computation is designed for the spiral CT images, where the Z-direction often has longer step length than those of X and Y-directions. If CT images are acquired using multi-slice CT with isotropic voxel size, the derivatives along diagonal directions in the X-Z and Y-Z planes must be also considered. For facilitating clarity, the illustrating example is provided in the Fig. 12.

The calculation of first-order derivative is described in the following equations. Let  $X_i$  be the intensity of *i*-th neighbor voxels and denote  $Y = \{y_i : i = 1, 2, ..., 5\}$  the vector of directional derivatives. Expression 1 shows how to calculate the Y. The neighborhood is defined in Fig. 13.

$$y_{1} = |x_{12} - 18x_{4} + 18x_{6} - x_{14}|/(12 \bullet \lambda_{xy}), \quad y_{2} = |x_{9} - 18x_{1} + 18x_{3} - x_{11}|/(12 \bullet \lambda_{y})$$

$$y_{3} = |x_{13} - 18x_{5} + 18x_{7} - x_{15}|/(12 \bullet \lambda_{xy}), \quad y_{4} = |x_{10} - 18x_{2} + 18x_{0} - x_{8}|/(12 \bullet \lambda_{x})$$

$$y_{5} = |x_{16} - x_{17}|/(2 \bullet \lambda_{z})$$

$$(1)$$

where  $\lambda$  is the step length along direction "\*." In the Equation 1, the central-5-point formula is employed to compute the first-order difference. However, the invention is not solely limited to it. Any other formula can be utilized for calculation of the first-order derivative in accounting of voxel size for any other applications.

The GFV is determined as  $min\{y_i, i = 1,2,3,4,5\}$ . The GFV is much less sensitive to the image noise than value of single directional derivative. Fig. 14 shows an example of GFV image for a virtual colonoscopy CT axial image.

By thresholding the GFV, the layer of boundary can be separated from the tissue homogeneous regions since the tissue homogeneous region has much lower GFV comparing to the tissue boundary where the intensity dramatically changes. The threshold can be pre-set for a specific application. As in the present invention for virtual colonoscopy CT images, the threshold can be set as 56. In other words, if the GFV of voxel is larger than 56, it is boundary voxels; Otherwise, the voxel is in a homogeneous tissue region. Usually, the nonuniform tagging in a region of tagged residue is slowly varying in intensity. That kind of intensity change does not contribute much to the value of GFV. That's why the GFV threshold cannot be too small. However, the threshold cannot be too large either since the extracted layer of boundary may not enclose the entire colon lumen. This is the result of non-uniform range of partial volume effect along different directions. If we reduce the threshold for GFV thresholding, we get better chance to enclose the colon lumen. However, that takes risk of over estimating the layer of boundary. To solve this problem, we fused the thresholded GFV result with the low-level voxel classification results. In other words, the information from both voxel classification and GFV are combined to achieve optimal estimation of layer of boundary. If voxel's GFV is larger than the preset threshold or the voxel is at the edge of region where all voxels are in the same cluster, the current voxel is a real boundary voxel.

The real boundary voxels form a boundary layer between air, colon wall, and tagged residue. The next step is to remove the tagged residue and to transform the boundary layer between colon wall and tagged residue into the boundary layer between air and colon wall. We call this procedure Boundary Layer Reconstruction. The reconstruction has to be implemented on each spatial separate region of tagged residue independently. In other words, the reconstruction method is the same, but the parameters of the method have to be adaptively adjusted to each region since each region of residue is in different tagging conditions.

For a selected region of tagged residue, the intensity of voxel that is not in the boundary layer will be set to the average air intensity of CT images, e.g. –850 HU. The average air intensity of CT image is denoted as T<sub>air</sub>. There are three kinds boundary voxels in the tagged residue region. The first is the air/residue boundary. The second is the wall/residue boundary. And the third is the boundary voxel that there are all three kinds of voxels around it. The first kind of boundary point is easy to detect. It is close to the air region and has greater GFV (e.g. GFV is larger than 500) since the intensity change from air to tagged region is the most dramatically in the CT image. When the region of the first kind of layer is determined, a conditional region growing algorithm will applied to the region. The region growing will expand the region into both tagged residue and air region. An example of the 6-connected 3-D region growing algorithm is described in the following:

```
Push the seed in the Queue

While( Queue is not empty )
{

Get the front voxel in the queue and denote its intensity as V<sub>c</sub>.

Label the current voxel as in the region.

Check the 6 closest 3-D neighbors of the current voxel.
{

If the neighbor is not labeled as in the region, denote the neighbor intensity as V<sub>n</sub>.

{

If( (V<sub>c</sub> > 150 HU) && (V<sub>c</sub>+50 HU ≤ V<sub>n</sub>) && ( V<sub>n</sub> < 350 HU )

{

Add this neighbor into the Queue

}

Else if( (V<sub>c</sub> < -100 HU) && (V<sub>c</sub>-50 HU ≥ V<sub>n</sub>) && (V<sub>n</sub> > T<sub>sir</sub>) )
{

Add this neighbor into the region

}

// End of checking this neighbor

} // End of checking all neighbors

pop the front voxel from the Queue.

} // While()
```

The expanded region of the first kind of layer will cover the area of partial volume layer right between the tagged residue and the air lumen. Intensities of voxels in the expanded region will be set to  $T_{sir}$  also. The expended region may also cover part of the layer of the third kind. The last step is to transform the layer of the  $2^{nd}$  and the  $3^{rd}$  into layer of air/tissue. Denote  $T_{ave}$  the average intensity of all voxels in the layer of the  $2^{nd}$  and the  $3^{rd}$ . Denote  $T_{tissue}$  the average intensity of soft tissue around colon lumen (e.g. -50 HU). Denote  $V_c$  the intensity of the current voxel in the

layer of  $2^{nd}$  and  $3^{rd}$  and  $V_f$  the transformed intensity value of the current voxel. The intensity transformation is defined by the following equation.

```
If ( V_c > T_{ave} ) {
V_t = T_{air};
}
Blue
V_f = \frac{T_{tissue} - T_{air}}{T_{ave} - T_{tissue}} * (T_{ave} - V_c) + T_{air}
}
```

After applying the transformation, the electronic cleansing on this region of tagged residue has been done. The Fig. 15 demonstrates a typical result of electronic cleansing. As shown in the picture, the thin layer of the attaching stool is also cleaned well.

We recruited 35 healthy male volunteers for evaluation of electronic cleansing method with laxative free bowel preparation. Among the volunteers, 22 are between 18 and 30 years old, 6 are between 31 and 40, and 7 are between 41 and 50. All volunteers are with no personal or family history of colon cancer or polyps and all of them have not any gastrointestinal symptom. All volunteers followed at least one of those 3 BPs and underwent virtual colonscopy CT scan. Some of them followed 2 or all 3 of those 3 BPs and underwent virtual colonscopy CT scan in different days. When the volunteers underwent CT scan, they got 2 series of scan at a time. One was for supine and another was for prone position. There were 60 cases, 120 series of scan, in total. An experienced radiologist reviewed electronic cleansing results of all cases utilizing the Viatronix V3D-Colon Module System. It was assumed that there was no polyp in all those healthy volunteers' colon. The radiologist reviewed electronic cleansed axial slices first to find any non-cleansed regions. Then, he flied through the 3-D colonic lumen model in order to search any cleansing artifact. There were two kinds of artifacts. One was error removal of part of haustral fold. Another was the residue material with irregular shape in the lumen. The latter was due to inadequate cleansing. Finally, he reviewed the original multi-planar reformatted images for imaging artifacts unrelated to electronic cleansing which might be related to colon motion or scanner protocol. The radiologist was allowed to correlate supine and prone images to confirm the cleansing error. A case was recorded as inadequate cleansing if one of the following situations occurred.

- 1. Inadequate cleansing situation 1: There was at least one tagged region larger than 5 mm in diameter not being cleansed. (This usually was found when reviewed the cleansed axial images).
- 2. Inadequate cleansing situation 2: There was at least one part of haustral fold thicker than 4 mm being removed. (This usually was found when flied through the 3-D lumen model).
- 3. Inadequate cleansing situation 3: There were two or more areas of residual irregular material with their shortest diameter larger than 4mm in the colonic lumen. (This usually were found when flied through the 3-D lumen model)

Fig. 16 lists the evaluation result. By applying Fisher test, it was concluded that BP 1 was significantly better than both BP 2 (p = 0.01) and BP 3 (p = 0.001). BP 2 had not significantly different from the BP 3. This result showed

that the BP 1, the fluid diet BP, provided promising results with adequate tagging quality and electronic cleansing result all the time. The BP 1 may be the most feasible laxative free bowel preparation for virtual colonsocpy.

The forgoing description focuses on oral intake of contrast agent. It is obvious that the electronic cleansing method is not limited to the oral intake contrast agent. It should work also for other kinds of contrast agent administration unless the tagged material has intensity difference to its surrounding tissues. The intravenous injection procedure for virtual clonoscopy is an example. The contrast agent should not limited to medicine or specific chemicals either. It could be normal food or even natural water. For example, as shown in the [3], the natural water was utilized as the contrast agent for imaging bladder. In that application, the organ preparation time is only half an hour. More general, the electronic cleansing methods can be applied to any modality of images has the same feature in intensity. The maximum directional gradient approach can also be applied to any modality of medical images for extraction of tissue boundary layer with partial volume effects.

The 3-D endoluminal view that was used in Viatronix V3D Colon Module was based on the volume rendering of the layer between the air and colon wall tissue. By reconstruction of the air/wall layer from the original residue/wall layer, it became feasible that the same rendering color map was able to apply to render entire cleansed colon lumen. The layer reconstruction was some kind of intensity transformation on the original volumetric data. This transformation was based on the information of both residue tagging and anatomy. The result of the cleansed volumetric images was the fusion of available information of both tagged residue and anatomy. In general, the fusion of information can be any kind of information that is related to the clinical application. The method of fusion of volumetric images with feature data is not limited to boundary layer reconstruction. The fusion method can be any transformation that maps feature data into the intensity of the volumetric images. The rendering method can also adaptively design for specific fusion technology and clinical requirement. The direct method is to adjust the color map of rendering to allow different intensity range showing different color in the endoluminal view.

In the application of laxative free virtual colonoscopy, computer assisted polyp detection (CAPD) results can be fused with the original volumetric images. The result of CAPD usually is a list of labeled regions for suspicious polyps and a list of likelihood that the region is related to be real polyp. One example of method to fuse the CAPD result with original volumetric images is to color code the suspicious regions in the volumetric images. Fig. 17 demonstrates an example. Another approach to fused the CAPD results is to transform the intensity of suspicious region into certain range and adjust the color map for volume rendering to allow the suspicious region show different color from surrounding normal tissue in the 3-D endoluminal view. For example, the clean colon wall and the wall coated with tagged stool can be shown in different color by adjusting the volume rendering color maps. For MRI virtual cystoscopy [3], the bladder wall region could be extracted with presented maximum gradient approach. Then, the texture analysis can be applied to the bladder wall region. The texture indexes associated to voxels can be mapped back to the original volumetric data and volume rendering in the endoluminal view for facilitating detection of abnormality in the bladder wall.

The foregoing merely illustrates the principles of the invention. It will thus be appreciated that those skilled in the art will be able to devise numerous systems, apparatus and methods which, although not explicitly shown or described herein, embody the principles of the invention and are thus within the spirit and scope of the invention as defined by its claims.

For example, the methods and systems described herein could be applied to virtually examine an animal or inanimate object. Besides the stated uses in the medical field, applications of the technique could be used to detect the contents of sealed objects which cannot be open or some of contents can be dyed to show contrast in the images.

#### What is claimed is:

1. A method for preparing volumetric images for virtual endoscopy was composite of the following steps:

Organ preparation;

Imaging the organ;

Processing acquired images;

Displaying images.

- 2. The method as recited in claim 1 wherein the organ preparation represented administration of contrast agents.
- 3. The method as recited in claim 1 wherein the organ preparation represented patient following specific diet.
- 4. The method as recited in claim 1 wherein the organ preparation represented combination of administrating contrast agent and patient following specific diet.
- 5. The method as recited in claim 2 wherein the contrast agent was taken orally.
- 6. The method as recited in claim 5 wherein the contrast agent was natural water that is taken by the patient half an hour prior to the MRI scan. The urine in the bladder forms a natural and non-invasive contrast agent.
- 7. The method as recited in claim 2 wherein the contrast agent was injected through intravenous.
- 8. The method as recited in claim 2 wherein the contrast agent was forced in the organ through a catherter.
- 9. The method as recited in claim 1 wherein the imaging of the organ was undergone with CT scanner.
- 10. The method as recited in claim 1 wherein the imaging of the organ was undergone with MRI.
- 11. The method as recited in claim 1 processing acquired images was composite of steps:

Segmentation of regions of interest;

Feature analysis to the regions of clinical interest;

Fusion of original volumetric images and feature data.

- 12. The method as recited in claim 11 wherein the segmentation of region of interest represented applying a low-level voxel classification first, then using anatomical knowledge to extract regions of clinical interest and to remove unrelated regions in the volumetric images.
- 13. The method as recited in claim 12 wherein the extraction of regions of clinical interest represented the boundary layer extraction that covers the wall of the hollow organ.

- 14. The method as recited in claim 12 wherein the extraction of region of clinical interest represented the extraction of the region of the lumen of the organ.
- 15. The method as recited in claim 11 wherein the feature analysis represented the computer assisted detection of abnormalities in the region of interest.
- 16. The method as recited in claim 11 wherein the fusion of original volumetric images and feature data represented the detection of the tagged region and replace the tagged region with other intensities of other tissue or material in the regions of interest.
- 17. The method as recited in claim 11 wherein the fusion of original volumetric images and feature data represented the reconstruction of extracted boundary layer into another kind of boundary layer in the volumetric images.
- 18. The method as recited in claim 11 wherein the fusion of original volumetric images and feature data represented transforming intensities of voxels in the original volumetric images that were in the detected suspicious abnormal regions into specific intensity range.
- 19. The method as recited in claim 1 wherein displaying images represented volume rendering the volumetric images directly.
- 20. The method as recited in claim 1 wherein displaying images represented volume rendering the fused volumetric images that was created as stated in claim 11.
- 21. The method as recited in claim 1 wherein displaying images represented volume rendering the inner surface of the hollow organ based on either the original volumetric images or the fused volumetric images.
- 22. The method as recited in claim 21 where volume rendering the inner surface of the organ represented adjusting the color map to display the suspicious area with different color within the context of the inner surface of the organ.
- 23. A method for laxative free virtual colonoscopy was composite of the following steps:

Bowel preparation without applying both oral laxative and suppository;

Distending colon with specific material;

Imaging colon;

Processing acquired images.

Displaying images.

- 24. The method as recited in claim 23 wherein the bowel preparation represented administrating contrast agent.
- 25. The method as recited in claim 24 wherein contrast agents were barium sulfate, diatrizoate meglumine, and diatrizoate sodium solution.
- 26. The method as recited in claim 23 wherein the bowel preparation represented the combination of administrating of contrast agent and patient following specific diet.
- 27. The method as recited in claim 26 wherein specific diet was any diet with low residue foods.
- 28. The method as recited in claim 23 wherein distending colon with room air.
- 29. The method as recited in claim 23 wherein distending colon with CO<sub>2</sub>.
- 30. The method as recited in claim 23 wherein imaging colon with CT or MRI.

- 31. The method as recited in claim 23 wherein processing acquired images represented electronic cleansing tagged residue in the colon lumen.
- 32. The method as recited in claim 23 wherein processing acquired images was composite of steps:

Electronic cleansing tagged residue in the colon lumen;

Automatic detection of suspicious polyps;

Fusing polyp detection results into the electronic cleansed volumetric images.

33. The method as recited in claim 31 wherein electronic cleansing tagged residue in the colon lumen was composite of steps:

Segmentation of region of interest;

Classification of tagged regions;

Extraction of boundary layers between tagged residue and other tissues or materials;

Cleansing tagged residue in the volumetric images;

Transforming the boundary layer between tagged residue and colon wall into the boundary layer between air and colon wall in the volumetric images.

- 34. The method as recited in claim 33 wherein segmentation of region of interest represented the detection and removal of regions of both lung and bone in the volumetric images by applying low-level voxel classification and using knowledge of anatomy.
- 35. The method as recited in claim 33 wherein segmentation of region of interest represented the extraction of both air and tagged residue filled colon lumen by applying low-level voxel classification and using knowledge of anatomy.
- 36. The method as recited in claim 33 wherein the classification of tagged regions represented the low-level voxel classification was applied to all tagged residue regions in the volumetric images. The number of created clusters was preset based on empirical knowledge.
- 37. The method as recited in claim 33 wherein the extraction of boundary layers between tagged residue and other tissues or materials represented applying maximum directional gradient approach to extract boundary layer.
- 38. The method as recited in claim 33 wherein the cleansing tagged residue in the volumetric images represented setting the intensity of non-boundary layer voxel in the tagged residue region into the average intensity of air voxels in the volumetric images.
- 39. The method as recited in claim 33 wherein the transformation from the residue/wall boundary layer into the air/wall boundary layer was a linear transformation.
- 40. The method as recited in claim 39 wherein the linear transforming was with a penalty term based on the maximum directional gradient value.
- 41. The method as recited in claim 23 wherein the displaying images represented volume rendering the inner surface of the colon lumen with specific color map that could show colon wall voxels in different color if they were in different intensity ranges.
- 42. A method for laxative free virtual colonscopy was composite of the following steps:

Administration of 4 doses of barium sulfates and 2 doses of Gastroview in 36 hours period prior to the CT scan;

Patient following the fluid diet (foods are listed in Fig. 6) in the last 24 hours prior to the scan and following a suggested food list (See Fig. 4) in 2 days before the last 24 hours;

Imaging the air/CO2 filled colon in both supine and prone position;

Segmentation of region of interest;

Electronic cleansing of tagged residue in the colon lumen;

Fusion the results of automatic polyp detection if it was available;

Display the fused volumetric images.

43. A method for laxative free virtual colonscopy was composite of the following steps:

Administration of 4 doses of barium sulfates and 2 doses of Gastroview in 36 hours period prior to the CT scan;

Patient taking foods from a given Meal Kit (foods are listed in Fig. 5) in the last 24 hours prior to the scan and following a suggested food list (See Fig. 4) in 2 days before the last 24 hours;

Imaging the air/CO2 filled colon in both supine and prone position;

Segmentation of region of interest;

Electronic cleansing of tagged residue in the colon lumen;

Fusion the results of automatic polyp detection if it was available;

Display the volumetric images.

44. A method for laxative free virtual colonscopy was composite of the following steps:

Administration of 4 doses of barium sulfates and 2 doses of Gastroview in 36 hours period prior to the CT scan;

Patient following a suggested food list (See Fig. 4) in 3 days prior to the CT scan;

Imaging the air/CO2 filled colon in both supine and prone position;

Segmentation of region of interest;

Electronic cleansing of tagged residue in the colon lumen;

Fusion the results of automatic polyp detection if it was available;

Display the volumetric images.

45. A method for electronic cleansing tagged residue in a hollow organ for virtual endoscopy was composite of the following steps:

Segmentation of the regions of tagged residue;

Extraction of the boundary layer between tagged residue and wall of the organ;

Transforming the intensity of voxels in the tagged region but not in the boundary layer into a preset intensity;

Transforming the boundary layer between tagged residue and wall of the organ into another known type of boundary layer in the volumetric images;

Displaying the transformed volumetric images.

- 46. A method as recited in claim 45 wherein Segmentation of the region of tagged residue represented a low-level voxel classification.
- 47. A method as recited in claim 45 wherein extraction of the boundary layer between tagged residue and wall of the organ represented the maximum directional gradient approach.
- 48. A method as recited in claim 47 wherein the maximum directional gradient approach was composite of steps:

Selecting several discrete directions in the image grid space to uniformly cover all directions;

Computing the first-order derivative along each selected direction for all voxels in the volumetric images;

Determining the maximum absolute value of all directional first-order derivative for all voxels in the volumetric images. This was the so-called maximum directional gradient value (MDGV).

Thresholding MDGVs. If the MDGV of a voxel was larger than a preset threshold, this voxel was in the boundary layer.

- 49. The method recited in claim 48 wherein the computing the first-order derivative along each selected direction for each voxel must consider the voxel sizes along different directions.
- 50. The method recited in claim 48 wherein the computing the first-order derivative along each selected direction for each voxel utilized different equation with longer range for reducing the affection of noise.
- 51. The method recited in claim 48 wherein the preset threshold for thresholding MDGV was determined by computing voxels in a typical tissue boundary region in the volumetric images.
- 52. The method recited in claim 45 wherein the transforming non-boundary tagged voxels into a preset intensity, the preset intensity was the average intensity of a specific tissue or material in the volumetric images.
- 53. The method recited in claim 45 wherein the transformation of the boundary tagged voxels was a linear transformation that flips intensity.
- 54. The method recited in claim 45 wherein displaying the transformed volumetric images represented the volume rendering the inner surface of the hollow organ. The viewing point was located in the cleansed lumen.
- 55. A method for extraction of layer of tissue or material boundary that have partial volume effect from medical images was composite of steps:

Segmentation of the regions of interest;

Extraction of the boundary layer;

Transform the non-boundary layer voxels into a preset intensity in the volumetric images;

Transform the boundary layer voxels into certain range of intensities in the volumetric images;

Displaying the transformed volumetric images.

- 56. The method recited in claim 54 wherein the extraction of the boundary layer represented the maximum directional gradient approach.
- 57. The method recited in claim 54 wherein the transformation of the boundary layer voxels represented replacing the intensity with the MDGV of the voxel.

- 58. The method recited in claim 54 wherein the displaying of the transformed volumetric images represented the volume rendering of the transformed volumetric images.
- 59. The method recited in claim 57 wherein the color map of volume rendering was adjusted to show different intensity range in different color.

# Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US04/033888

International filing date: 12 October 2004 (12.10.2004)

Document type: Certified copy of priority document

Document details: Country/Office: US

Number: 60/510,237

Filing date: 10 October 2003 (10.10.2003)

Date of receipt at the International Bureau: 12 November 2004 (12.11.2004)

Remark: Priority document submitted or transmitted to the International Bureau in

compliance with Rule 17.1(a) or (b)



## This Page is Inserted by IFW Indexing and Scanning Operations and is not part of the Official Record

#### **BEST AVAILABLE IMAGES**

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

☐ BLACK BORDERS
$\square$ image cut off at top, bottom or sides
FADED TEXT OR DRAWING
☑ BLURRED OR ILLEGIBLE TEXT OR DRAWING
☐ SKEWED/SLANTED IMAGES
☐ COLOR OR BLACK AND WHITE PHOTOGRAPHS
☐ GRAY SCALE DOCUMENTS
☐ LINES OR MARKS ON ORIGINAL DOCUMENT
☐ REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY

#### IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.